



---

# Bencox<sup>®</sup>

## Mirabo Z Cup


---



# INNOVATIVE

- Secure Locking Mechanism
- Wear Mechanism
- Laser 3D Coating Technology






# » Secure Locking Mechanism


## Locking with PE Liner



The patented\* locking mechanism allows for secure fixation of the cup and liner with excellent resistance to push-out, rotation and tilting.



12 anti-rotational tabs  
\* Patent no. : 9,144,497 (U.S.A)  
2 674 132 (EUROPE)  
ZL201180067109.1 (CHINA)  
10-1268106 (KOREA)  
10-1279095 (KOREA)



# » Secure Locking Mechanism


## Locking with Ceramic Liner

---

Cup and ceramic insert provide a secure taper locking mechanism while providing a self-locating feature to ease the insert placement.

Increased taper length prevents mal-positioning of the ceramic liner.


• Taper Lock



# » Surface Treatment

## TPS Coating


The microporous pure titanium coating with pores of 50-200 $\mu$ m diameter and more than 30% porosity leads to direct bone apposition for secondary stability.



Mag. 100x

## Laser 3D Coating

Laser 3D coating technique produces layers of titanium powder using laser beam. The coating with pores of 200-500 $\mu$ m diameter and 50~70% porosity leads to bone-in-growth. Increased roughness from the Laser 3D coating allows about 1.4 times stronger initial fixation than Ti Plasma Spray coating.



Mag. 100x




## » Secure Peripheral Fixation of Cup


Slightly flattened dome establishes physiological load transfer throughout the periphery of the cup.



FE analysis model




Secure Initial Fixation



## Enhanced Wear Characteristics

In-vitro tests showed that Mirabo PE liner exhibited 99.8% less wear than the conventional PE liners.



Standard type



15° Elevated type


# » Bencox<sup>®</sup> Mirabo PE Liner

Bencox<sup>®</sup>  
PE Liner


## Highly Crosslinked Polyethylene Manufacturing Process

- Compression molded UHMWPE (GUR 1050) rod stock material.
- Gamma irradiated at 10Mrad at room temperature to facilitate cross-linking.
- Following irradiation, the rods are heated above the melting point to eliminate residual free radicals and form additional cross-links.
- PE liners are machined, cleaned, packaged, and EtO sterilized.

10Mrad  
Gamma  
Irradiation




UHMWPE  
(GUR 1050)



Cross linking with  
free radicals present


Remelting



Free radicals  
eliminated

Ordering Information

| Part No.               |          |                               |          | Dimensions |            |            |            |           |
|------------------------|----------|-------------------------------|----------|------------|------------|------------|------------|-----------|
| Subject Device         |          | Predicate Device<br>[K172806] |          | (D1)<br>mm | (D2)<br>mm | (H1)<br>mm | (H2)<br>mm | (T)<br>mm |
| BENCOX<br>Mirabo Z Cup | Hole Qty | BENCOX<br>Mirabo Cup          | Hole Qty |            |            |            |            |           |
| H1.A23.1044            | 3        | H1.A11.1044                   | 3        | 44         | 35         | 23.5       | 19.1       | 3.2±0.5   |
| H1.A23.1046            | 3        | H1.A11.1046                   | 3        | 46         | 37         | 24.8       | 19.8       | 3.6±0.5   |
| H1.A23.1048            | 3        | H1.A11.1048                   | 3        | 48         | 39         | 24.8       | 19.8       | 3.2±0.5   |
| H1.A23.1050            | 3        | H1.A11.1050                   | 3        | 50         | 39         | 25.8       | 19.8       | 4.2±0.5   |
| H1.A23.1052            | 3        | H1.A11.1052                   | 3        | 52         | 44         | 27.9       | 22.3       | 3.5±0.5   |
| H1.A23.1054            | 3        | H1.A11.1054                   | 3        | 54         | 44         | 27.6       | 22.3       | 3.4±0.5   |
| H1.A23.1056            | 3        | H1.A11.1056                   | 3        | 56         | 44         | 28.0       | 22.3       | 3.9±0.5   |
| H1.A23.2056            | 3        | H1.A11.2056                   | 3        | 56         | 48         | 28.5       | 23.4       | 3.6±0.5   |
| H1.A23.1058            | 3        | H1.A11.1058                   | 3        | 58         | 48         | 28.9       | 23.4       | 3.6±0.5   |
| H1.A23.1060            | 3        | H1.A11.1060                   | 3        | 60         | 48         | 29.9       | 23.4       | 4.6±0.5   |
| H1.A23.1062            | 3        | H1.A11.1062                   | 3        | 62         | 48         | 30.8       | 23.4       | 5.6±0.5   |
| H1.A23.1064            | 3        | H1.A11.1064                   | 3        | 64         | 52         | 31.8       | 25.3       | 4.9±0.5   |
| H1.A23.1066            | 3        | H1.A11.1066                   | 3        | 66         | 52         | 32.8       | 25.3       | 5.9±0.5   |
| H1.A23.1068            | 3        | H1.A11.1068                   | 3        | 68         | 52         | 33.7       | 25.3       | 6.9±0.5   |





Bencox<sup>®</sup>

Mirabo Z Cup



|                        |                                                                                                   |
|------------------------|---------------------------------------------------------------------------------------------------|
| Manufacturing Plant    | 12, Yeongsanhong 1-gil, Ipjang-myeon, Seobuk-gu, Cheonan-si, Chungcheongnam-do, Republic of Korea |
| Seoul Office           | 33-2, Banpo-daero 20-gil, Seocho-gu, Seoul, Republic of Korea                                     |
| Corentec America, Inc. | 23181 Verdugo Dr. #104A, Laguna Hills, CA 92653                                                   |

[www.corentec.com](http://www.corentec.com)  
CRSC-0107-M01(KR00)