


OP set variable angle-stable proximal humerus plate

8 head holes, for screws 3.5/4.0 mm, small fragment, titanium, complete set item no. 19.3673.00



item no. 5.790.84



OP set complete item no. 19.3673.00

item description	item no.	quantity
Implants		
<u>Variable angle-stable humerus plate with combination holes in shank</u>		
8/4 holes, for screws Ø 4.0/3.5; length 93 mm	5.790.84	1
8/6 holes, for screws Ø 4.0/3.5; length 125 mm	5.790.86	1
8/8 holes, for screws Ø 4.0/3.5; length 156 mm	5.790.88	1
plates available on request up to 14 holes (length 253 mm)		

Instruments		
drill bit for quick coupling, Ø 2.5 x 155, spiral length 60, double spiral	2.904.22	1
drill guide for angle stable screwing, length 68, small fragment	2.977.04	2
hex screwdriver w. handle f. screws Ø 2.7 - Ø 4.0, AF2.5, conical	2.9406.25	1
gauge w. clasp f. screws w. con. head thread Ø 3.5, Ø 4.0, MR 70	2.953.70	1
guide block f. prox. SF humerus plates, 8 holes in head, flat, PEEK	10.179.10	1
sleeve for guide block for proximal small fragment humerus plates	10.173.01	2
screw forceps, self holding	2.954.01	1
drill guide w. handle f. var. angle-stable screwing, SF, length 68	2.977.09	1
screw driver for guide block	10.173.05	1



10.179.10



10.173.01



10.173.05

Container		
perforated autoclavable container w. inset f. instruments for humerus-set	19.366.00	1

Wire		
Kirschner wire with trocar point and round end, Ø1.8 x 150, st.st.	6.031.18	5



- Variable angle-stable holes for setting angle-stable screws
- Wire wrap holes
- K-Wire holes
- Combination holes (can be set with standard or angle-stable screws)
- Longhole (for positioning of the plate with standard screws)
- holes for guide block

item description	item no.	quantity
angle-stable cortical screws		
conical head thread, fully threaded, Ø3.5 mm, self-tapping, titanium		
length 22 mm	3.133.22	3
length 24 mm	3.133.24	5
length 26 mm	3.133.26	5
length 28 mm	3.133.28	5
length 30 mm	3.133.30	5
length 32 mm	3.133.32	5
length 34 mm	3.133.34	5
length 36 mm	3.133.36	5
length 38 mm	3.133.38	5
length 40 mm	3.133.40	5
length 42 mm	3.133.42	5



cortical screws		
fully threaded, self-tapping, titanium Ø 3,5 mm		
length 22 mm	3.132.22	5
length 24 mm	3.132.24	5
length 26 mm	3.132.26	5
length 28 mm	3.132.28	5
length 30 mm	3.132.30	5
length 32 mm	3.132.32	5
length 34 mm	3.132.34	5
length 36 mm	3.132.36	5
length 38 mm	3.132.38	5
length 40 mm	3.132.40	5
length 42 mm	3.132.42	5



angle-stable cancellous screws		
conical head thread, fully threaded, titanium self-tapping, Ø 4.0 mm		
length 34 mm	3.147.34	5
length 36 mm	3.147.36	5
length 38 mm	3.147.38	5
length 40 mm	3.147.40	5
length 42 mm	3.147.42	5
length 44 mm	3.147.44	5
length 46 mm	3.147.46	5
length 48 mm	3.147.48	5
length 50 mm	3.147.50	5
length 52 mm	3.147.52	5
length 54 mm	3.147.54	5
length 56 mm	3.147.56	5



OP set variable angle-stable proximal humerus plate

8 head holes, for screws 3.5/4.0 mm, small fragment, titanium, complete set item no. 19.3673.00

